

Early Therapeutic Discovery in Oncology and Beyond: Challenges and Opportunities

Anton Simeonov, Ph.D.


Scientific Director, National Center for Advancing Translational Sciences (NCATS), National Institutes of Health (NIH)

TRACO Lecture
September 28, 2020



Drug Discovery Process

Steps in the Drug Discovery Process




Range of screening assays

Range of Screening Assays

Extent of reductionism →


Phenotype
*(Image-based
HCS, GFP, etc)*



Pathway
*(Reporters, e.g., luciferase, β -
lactamase)*




Protein
(Enzyme readouts, interactions, etc)



High throughput screening


Robotics & Informatics




Microliter Dispensing



Laser Cytometry



Microscopy




HTS is a standard step in the drug discovery process but has remained problem-ridden.



National Center
for Advancing
Translational Sciences

Cell based HTS assays

149 Cancer Relevant Cell-Based HTS Assays from PubChem



Coussons, N. P., Braisted, J. C., Peryea, T., Sittampalam, S. G., Simionov, A. and Hall, M. D. **Small Molecule Screens: A Gateway to Cancer Therapeutic Agents with Case Studies of FDA-Approved Drugs** *Pharmacological Reviews*, October 2017, 69 (4) 479-496



Assay choice

Important Considerations for Choosing an Assay

- **Assay expense**
 - Cost per well
 - Disposal cost(s)

Important considerations

Important Considerations for Choosing an Assay

- **Assay expense**
 - Cost per well
 - Disposal cost(s)
- **Available instrumentation**
 - Select the best possible assays based on the available instrumentation

Assay throughput



Important Considerations for Choosing an Assay

- **Assay expense**
 - Cost per well
 - Disposal cost(s)
- **Available instrumentation**
 - Select the best possible assays based on the available instrumentation
- **Assay throughput**
 - Miniaturization reduces the cost per well



Assay miniaturization

Assay Miniaturization Saves Time and Reagents



	96	384	1536
Plates per 100,000 compounds:	1,042	261	66
Assay volume (uL):	50-200	30-50	2-8
Adherent cell seeding density:	~10,000	~2,000	~500

Jordan, Shane R. 'Continuous High Content Phenotypic Screening: Success Factors in Drug Discovery.' InTech, 2018.

Important considerations

Important Considerations for Choosing an Assay

- **Assay expense**
 - Cost per well
 - Disposal cost(s)
- **Available instrumentation**
 - Select the best possible assays based on the available instrumentation
- **Assay throughput**
 - Miniaturization reduces the cost per well
- **Ability to multiplex**
 - Can the response be measured by a single parameter; is multiparametric output possible?
 - Increased data per sample
 - Can guide hit selection by differentiating selectivity among related targets
 - Can distinguish pathway inhibition from cytotoxicity in a cell-based assay



Important considerations

Important Considerations for Choosing an Assay

- **Assay expense**
 - Cost per well
 - Disposal cost(s)
- **Available instrumentation**
 - Select the best possible assays based on the available instrumentation
- **Assay throughput**
 - Miniaturization reduces the cost per well
- **Ability to multiplex**
 - Can the response be measured by a single parameter; is multiparametric output possible?
 - Increased data per sample
 - Can guide hit selection by differentiating selectivity among related targets
 - Can distinguish pathway inhibition from cytotoxicity in a cell-based assay
- **Reagents**
 - Stability for hours is important
 - Consistency is critical (ideally obtain a large quantity from a single lot)
 - All reagents need to be validated (cell lines, antibodies, enzymatic purity, etc.)



Important considerations


Important Considerations for Choosing an Assay

- **Assay expense**
 - Cost per well
 - Disposal cost(s)
- **Available instrumentation**
 - Select the best possible assays based on the available instrumentation
- **Assay throughput**
 - Miniaturization reduces the cost per well
- **Ability to multiplex**
 - Can the response be measured by a single parameter; is multiparametric output possible?
 - Increased data per sample
 - Can guide hit selection by differentiating selectivity among related targets
 - Can distinguish pathway inhibition from cytotoxicity in a cell-based assay
- **Reagents**
 - Stability for hours is important
 - Consistency is critical (ideally obtain a large quantity from a single lot)
 - All reagents need to be validated (cell lines, antibodies, enzymatic purity, etc.)
- **Potential for assay interference**
 - Fluorescent compounds can interfere with fluorescent readouts
 - Colored compounds might interfere with luminescence



Fluorescence spectroscopic profiling



Fluorescence Spectroscopic Profiling of Compound Libraries



Simionov, A., Jadhav, A., Thomas, C.J., Wang, Y., Huang, R., Southall, N.T., Shinn, P., Smith, J., Austin, C.P., Auld, D.S. and Ingleso, J., 2008. Fluorescence spectroscopic profiling of compound libraries. *Journal of Medicinal Chemistry*, 51(8), 2363-2371.

Assay tolerance

Determination of Assay Tolerance to DMSO/Vehicle is Important



Yasgar A., Jadhav A., Simoonov A., Coussons N.P., **AlphaScreen-Based Assays: Ultra-High-Throughput Screening for Small Molecule Inhibitors of Challenging Enzymes and Protein-Protein Interactions**. *Methods Mol Biol.* 2016;1439:77-98.



Important considerations

Important Considerations for Choosing an Assay

- Homogenous assay format is preferred for screening
 - Add reagents, mix and measure (no solution removal or wash steps)
 - Automation friendly
 - Reduces variability
 - Decreases hands-on time
 - Improves reproducibility

Important considerations

Important Considerations for Choosing an Assay

- Homogenous assay format is preferred for screening
 - Add reagents, mix and measure (no solution removal or wash steps)
 - Automation friendly
 - Reduces variability
 - Decreases hands-on time
 - Improves reproducibility
- Time required for assay
 - Off-line reagent preparation
 - Is temperature equilibration required
 - Actual assay time
 - Kinetic versus end point read
 - Time required for data analysis and record keeping



Important considerations

Important Considerations for Choosing an Assay

- Homogenous assay format is preferred for screening
 - Add reagents, mix and measure (no solution removal or wash steps)
 - Automation friendly
 - Reduces variability
 - Decreases hands-on time
 - Improves reproducibility
- Time required for assay
 - Off-line reagent preparation
 - Is temperature equilibration required
 - Actual assay time
 - Kinetic versus end point read
 - Time required for data analysis and record keeping
- Signal stability
 - Does the response occur rapidly or within a few minutes or hours?
 - Longer signal stability allows for flexibility in automated systems
 - Longer signal stability minimizes differences among plates within a stack



Important considerations


Important Considerations for Choosing an Assay

- Homogenous assay format is preferred for screening
 - Add reagents, mix and measure (no solution removal or wash steps)
 - Automation friendly
 - Reduces variability
 - Decreases hands-on time
 - Improves reproducibility
- Time required for assay
 - Off-line reagent preparation
 - Is temperature equilibration required
 - Actual assay time
 - Kinetic versus end point read
 - Time required for data analysis and record keeping
- Signal stability
 - Does the response occur rapidly or within a few minutes or hours?
 - Longer signal stability allows for flexibility in automated systems
 - Longer signal stability minimizes differences among stacks of plates
- **Assay Sensitivity**
 - Choice of readouts is important
 - Colorimetric < fluorescent < luminescent



Assay suitability

Evaluating Assay Suitability for Screening




$S/N = 12$

$S/B = 5$

$z = 0.5$

$$S/N = \frac{\text{mean signal} - \text{mean background}}{\text{standard deviation of background}}$$

$$S/B = \frac{\text{mean signal}}{\text{mean background}}$$



$S/N = 27$

$S/B = 10$

$z = 0.1$


$$Z = 1 - \frac{3SD \text{ of sample} + 3SD \text{ of control}}{|\text{mean of sample} - \text{mean of control}|}$$

A Simple Statistical Parameter for Use in Evaluation and Validation of High Throughput Screening Assays. Zhang JH, Chung TD, Oldenburg KR. J Biomol Screen. 1999;4(2):67-73.



Assay suitability

Evaluating Assay Suitability for Screening



$$Z = 1 - \frac{3SD \text{ of sample} + 3SD \text{ of control}}{|\text{mean of sample} - \text{mean of control}|}$$


Z-factor value	Structure of assay	Related to screening
1	SD = 0 (no variation), or the dynamic range $\rightarrow \infty$	An ideal assay
$1 > Z \geq 0.5$	Separation band is large	An excellent assay
$0.5 > Z > 0$	Separation band is small	A double assay
0	No separation band, the sample signal variation and control signal variation bands touch	A "yes/no" type assay
<0	No separation band, the sample signal variation and control signal variation bands overlap	Screening essentially impossible

A. B. and S. A. S. C. Parameter for Use in Evaluation and Validation of High Throughput Screening Assays, Zhang L. L. Chung T. D. Odellburg K. L. C. Data Screen. 1999;4(2):67-70.



Improving early discovery

Improving the Process of Early Discovery: Quantitative High-Throughput Screening (qHTS)



- Conventional screening done at one concentration
 - Not appropriate for potency testing – “dose makes the poison”
- qHTS tests compounds assayed at **multiple** concentrations (range: 4 logs)
- Enabled by miniaturized assay volumes (2-8 μL per test) and informatics pipeline
- Generates *pharmacological actives* instead of statistical “hits”
 - Dramatically increases reliability
 - Dramatically reduces false positives and false negatives
- *To date, several hundred million datapoints from several hundred screens have been generated and deposited in the public domain.*

Medicinal chemistry

Medicinal Chemistry, an Integrated Process



Tier 1: Synthesis & validation



- Medicinal chemistry
- Purification
- In vitro ADME

Tier 2: Compound profile expansion

- Met ID/ CYP studies
- In vitro toxicology

Tier 3: Advanced Preclinical studies

- Formulation
- Scale-up
- In vivo PK/PD
- Preclinical toxicology



Small molecular probes

Example: small molecule probes to study cancer metabolism...with additional applications

- Tumors exhibit unique nutritional requirements. For example, pyruvate kinase M2 (PKM2) has been implicated in the Warburg Effect.
- A screen of PKM2, followed by medicinal chemistry optimization, produced several distinct molecules able to significantly activate the enzyme.
- These tool molecules have been used by numerous labs worldwide to study cancer biology and a subset are being further developed as therapies.

PKM2 activation may also be a way to modulate metabolism in pancreatic beta cells to manage diabetes:

ARTICLES

nature
medicine

Pyruvate kinase M2 activation may protect against the progression of diabetic glomerular pathology and mitochondrial dysfunction

Inhibition of Pyruvate Kinase M2 by Reactive Oxygen Species Contributes to Cellular Antioxidant Responses

Dimitrios Anastasiou,^{1,2} George Pasadounakis,^{1,2} John B. Baena,^{1,2} Matthew B. Rosen,⁴ Ben-Kang Jiang,⁴ Min Shen,⁴ Gary Bollinger,^{1,2} Aron T. Sasaki,^{1,2} Jesse M. Lincenko,^{1,2} Douglas S. Auld,^{1,2} Craig J. Thomas,² Matthew G. Vander Heiden,^{1,2} Lewis C. Cantley^{1,2,3}



Inhibition of Pyruvate Kinase M2 by Reactive Oxygen Species Contributes to Cellular Antioxidant Responses
Dimitrios Anastasiou et al.
Science 334, 1278 (2011);
DOI: 10.1126/science.1211485


Because IPMs are sufficiently general and because density dependence and environmental variation affect most populations, these conclusions are likely to extend to other systems. The construction and analysis of IPMs across a range of systems may provide support for this proposition. In addition to providing a tool to explore eco-evolutionary dynamics, IPMs have also been extended to include spatial variation and to identify evolutionarily stable strategies (21, 22), giving them potential to unify several subdisciplines of population biology, including population ecology, quantitative genetics, population genetics, and



NIH
National Center
for Advancing
Translational Sciences

Drug combinations



Translation Challenge: Rapid Discovery of Drug Combinations



National Center
for Advancing
Translational Sciences

Resistance

Application of Drug Combinations to Address Resistance



Drug resistance



Dissemination of technology: combination screening to overcome drug resistance in cancer cells

- ABC subtype of Diffuse Large B-Cell Lymphoma (ABC-DLBCL) has poor prognosis and response to treatment
- Ibrutinib is a BTK inhibitor that has activity against ABC DLBCL
- Study evaluated 459 drugs in combination with Ibrutinib
 - » 6 x 6 concentration-response “matrix blocks”, validation in 10 x 10 concentration-response matrix blocks
- DNA-damaging agents identified as synergizing with Ibrutinib in killing ABC DLBCL cell lines
- **Dissemination:**
 - » Protocols
 - » Source code for dispense

High-throughput combinatorial screening identifies drugs that cooperate with ibrutinib to kill activated B-cell-like diffuse large B-cell lymphoma cells

Lesley A. Mathews-Gibson^{1,2}, Rajarshi Guha^{1,2}, Paul Shien^{1,2}, Ryan M. Young^{1,2}, Jonathan M. Keller¹, Songbo Liu¹, Ian S. Goldfarb¹, Adam Yeager¹, Crystal McKinstry¹, Matthew B. Bower¹, Damien Y. Dumesnil¹, Jian-Kang Jiang¹, Sam Michael¹, Tim Mierzwski¹, Wenwei Huang¹, Martin J. Walsh¹, Bryan T. Mott¹, Parasima Patel^{1,2}, William Golder¹, David A. Winkler¹, Christopher A. Loda¹, Giuseppe Ruff¹, Ajit Keshav¹, Brian D. Pepper¹, Christopher R. Austin¹, Scott E. Marder¹, Anton Simakov¹, Marc Ferrer¹, Leslie M. Stauch^{1,2}, and Craig J. Thomas^{1,2}


¹Division of Precision Innovation, National Institutes of Health Chemical Genomics Center, National Center for Advancing Translational Sciences, Biomedical Research Center for Cancer Research, and ²Translational Therapeutics Program, Division of Cancer Treatment and Diagnosis, National Cancer Institute, National Institutes of Health, Bethesda, MD 20892 and ³Translational Therapeutics Program, Division of Cancer Treatment and Diagnosis, National Cancer Institute, National Institutes of Health, Bethesda, MD 20892



Rare brain tumors



Example: combination screening to tackle untreatable rare brain tumors

Panobinostat
- vs. -
MPE 5.0 library
(2450 agents)
Marizomib
- vs. -
MPE 5.0 library
(2450 agents)




Lin *et al.*, *Sci. Transl. Med.* **11**, eaaw0064 (2019) 20 November 2019

Collaboration between NCATS (Craig Thomas) and Stanford University (Michelle Monje)



Panobinostat
HDAC inh.
Approved
MPO score: 4.82
Limited CNS pen.

+



Marizomib
Proteasome inh.
Phase 1
MPO score: 5.5
CNS penetrant

The team is currently exploring a clinical trial of panobinostat together with marizomib in DMGs.

SCIENCE TRANSLATIONAL MEDICINE | RESEARCH ARTICLE


CANCER

Therapeutic strategies for diffuse midline glioma from high-throughput combination drug screening


3D models

Increasing the predictivity of *in vitro* assays: a continuum of 3D models of healthy and diseased tissues


2D




Spheroids




Organoids



Printed Tissues





Organ-on-a-chip



HTS compatibility


Physiological complexity



National Center
for Advancing
Translational Sciences


Tissue bioprinting

3D Tissue Bioprinting



Gel

+



Cells




Syringe



Printer

Hydrogel polymer is mixed with cells and loaded into syringe.



The printer "3D prints" the cell/gel mixture in a layer by layer approach.



Printed construct



1 day



1 week



2 weeks

The printed construct is incubated to allow the cells to form a tissue, and to enable proper cell differentiation.

Epidermis Functional activity analyses. Stem cell technologies

Layers of the Epidermis: native skin *versus* 3D-bioprinted skin

Native Skin



3D-Bioprinted Skin


Stratum corneum

**Stratum
granulosum**

Stratum spinosum

**Stratum
germinativum**

**Dermis:
Collagen and elastic
fibers**



<http://www.sjumed.edu/~eking2/intro/IN005b.htm>


3D Tissue model

www.oncotarget.com

Oncotarget, 2020, Vol. 11, (No. 27), pp: 2587-2596

Research Paper

A 3D biofabricated cutaneous squamous cell carcinoma tissue model with multi-channel confocal microscopy imaging biomarkers to quantify antitumor effects of chemotherapeutics in tissue




Collaboration between NCATS (Marc Ferrer) and Rockefeller University (Daniel Gareau)



Stem cell technologies


Enabling Advanced 3D Models and Regenerative Medicine through Stem Cell Technologies

NCATS Stem Cell Translation Laboratory:




Access to relevant human cell types

Sensory neurons (nociceptors) and other neuronal subtypes
Hepatocytes, etc




Advanced imaging technologies for functional cell characterization

High content confocal calcium imaging optogenetics




High-throughput electrophysiology methods

High density multi-electrode arrays
28400 electrode well




Measurement of signaling pathways, metabolism & specific targets

Cyclic AMP, PKA activity, CREB phosphorylation
energy metabolism



Longitudinal tracking of cell behavior

Multiple measurements over days, weeks or months



Combined single-cell transcriptomic & proteomic analyses

Drug response in individual cells



Assay development

**Where do I go for more
information about assay
development?**



Assay guidance manual

Sharing internal know-how: Assay Guidance Manual (47 chapters/ 1,338 printed pages)




Table of Contents

Preface	
Considerations for Early Phase Drug Discovery	1 Chapter
In Vitro Biochemical Assays	10 Chapters
In Vitro Cell-Based Assays	10 Chapters
In Vivo Assay Guidelines	2 Chapters
Assay Artifacts and Interferences	4 Chapters
Assay Validation, Operations and Quality Control	5 Chapters
Assay Technologies	2 Chapters
Instrumentation	2 Chapters
Pharmacokinetics and Drug Metabolism	1 Chapter
Glossary of Quantitative Biology Terms	1 Chapter

<https://ncats.nih.gov/agm-video>

August 2nd Videos

1. Assav, CP: Welcome to the Assay Guidance Manual (AGM) Workshop
2. Carosons, NP: Strategies for Assay Selection & Robust Biochemical Assays
3. Ass, T: Treating Cells as Reagents to Design Reproducible Screening Assays
4. Fink, CL: Assay Development Considerations for High Content Imaging
5. Avila, DS: Studies in Interference and Interference in Assay Interference
6. Dabke, J: Assay Interference by Chemical Reactivity
7. Chung, TD: Basic Assay Statistics, Data Analysis & Rules of Thumb
8. Damsky, M: Reproducibility & Differentiability of Potency Results
9. Swamishan, GS: Avoiding Artifacts & Interference in Assay Operations

March 20-21st Videos

1. Assav, CP: Welcome to the Assay Guidance Manual (AGM) Workshop
2. Carosons, NP: Robust Assays Define Success in Preclinical Research
3. Huang, M: Target Identification & Validation in Translational Discovery
4. Foley, T: Development & Validation of Cell-Based and Biochemical Assays
5. Ass, T: Treating Cells as Reagents to Design Reproducible Screening Assays
6. Fink, CL: Assay Development for HCS & Best Practices for 3D HCS
7. Kach, C: Mass Spectrometry for Drug Screening and Lead Optimization
8. Dabke, J: Bioassay Interference by Aggregation and Chemical Reactivity
9. Fink, CL: Lead Selection and Optimization by Medicinal Chemistry
10. Fink, CL: In Vitro Technological Testing Using a qHTS Platform
11. Ass, T: In Vitro Assessment of ADME Properties of Lead Compounds
12. Calk, SD: Statistical Design of Experiments for Assay Development
13. Gupta, R: Photos Application to Target Evaluation and Drug Discovery
14. Wodner, JR: Assay Operations: Keeping Assays Robust and Reproducible

Website: <https://ncats.nih.gov/exportofpreclinicalassays>

Email us: NCATS_AGM_Editorial@mail.nih.gov



Facebook: www.facebook.com/assayguide



LinkedIn: www.linkedin.com/company/ncats/7497344




National Center
for Advancing
Translational Sciences

Response to Covid-19

Response to COVID-19: OpenData Portal enables data and protocol sharing in near-real time

U.S. Department of Health and Human Services > National Institutes of Health > National Center for Advancing Translational Sciences



OpenData Portal

Home

OpenData Browser

Assays

Animal Models

Omics Efforts

Highlights

Resources ▾

OpenData | COVID-19

NCATS is generating a collection of datasets by screening a panel of SARS-CoV-2-related assays against all approved drugs.

These datasets, as well as the assay protocols used to generate them, are being made immediately available to the scientific community on this site as these screens are completed.

<https://opendata.ncats.nih.gov/covid19/>

Open data portal



bioRxiv

THE PREPRINT SERVER FOR BIOLOGY

An OpenData portal to share COVID-19 drug repurposing data in real time

Kyle R. Brimacombe, Tongan Zhao,  Richard T. Eastman, Xin Hu, Ke Wang, Mark Backus, Bolormaa Baljinnyam, Catherine Z. Chen, Lu Chen, Tara Eicher, Marc Ferrer, Ying Fu, Kirill Gorshkov, Hui Guo, Quinlin M. Hanson, Zina Itkin, Stephen C. Kales, Carleen Klumpp-Thomas, Emily M. Lee, Sam Michael, Tim Mierzwa, Andrew Patt, Manisha Pradhan, Alex Renn, Paul Shinn, Jonathan H. Shrimp, Amit Viraktamath, Kelli M. Wilson, Miao Xu, Alexey V. Zakharov, Wei Zhu, Wei Zheng, Anton Simeonov, Ewy A. Mathé, Donald C. Lo,  Matthew D. Hall, Min Shen

doi: <https://doi.org/10.1101/2020.06.04.135046>

This article is a preprint and has not been certified by peer review [what does this mean?].

Abstract

Full Text

Info/History

Metrics

 Preview PDF

Abstract

The National Center for Advancing Translational Sciences (NCATS) has developed an online open science data portal for its COVID-19 drug repurposing campaign – named OpenData – with the goal of making data across a range of SARS-CoV-2 related assays available in real-time. The

Multiple publicity, incl. featured on Francis Collins' blog post



Learn More About NCATS



Website: www.ncats.nih.gov



Facebook: facebook.com/ncats.nih.gov



Twitter: twitter.com/ncats_nih_gov



YouTube: youtube.com/user/ncatsmedia



E-Newsletter: ncats.nih.gov/news-and-events/e-news/e-news.html

Email us: info@ncats.nih.gov



National Center
for Advancing
Translational Sciences